

RNAscope HiPlex Experiment Checklist for Fixed FROZEN Samples				
Method: RNAscope HiPlex				
Experiment Name			Sample Info	
DATE:				
Sample Preparation Steps				
Prechill 4% PFA on ice or in 4 °C fridge				
1	OCT Removal	1X PBS	5min	RT
2	Drying	NA	30min	60°C
3	Post-fixation on ice	4% PFA, cold	15min	4°C
4	Dehydration			
		a 50% EtOH (50 ml)	5min	RT
		b 70% EtOH (50 ml)	5min	RT
		c 100% EtOH (50 ml)	5min	RT
		d 100% EtOH (50 ml)	5min	RT
		e 100% EtOH (50 ml)	5min	RT
Slides may be stored in 100% Ethanol at -20°C for 1 week				
Target Retrieval Using Temperature Controlled Heater				
Prewarm the Target Retrieval Buffer to >= 99°C, but do not boil for > 15 min before use				
5	Acclimate to 99°C	DW preheated in steamer	Briefly	99°C
6	Target Retrieval	Retrieval Buffer (kit)	5 min	99°C
7	Rinse	DW	15 sec	RT
8	Dehydration	100% EtOH	3 min	RT
9	Drying		5 min	60°C
10	Create hydrophobic barrier, let barrier dry completely	NA	5min	RT
Place a humidifying paper in humidity control tray and wet completely with distilled water				
11	Load slides into ACD EZ-Batch Slide Holder	NA	Briefly	RT
12	Protease Digestion	Protease III (kit)	15-30 min	40°C
13	Rinse X2	PBS	2 min	RT
HiPlex Protocol				
Round 1	Prepare 3 L 1 X Wash Buffer, 4 X SSC, 1 L 0.5% PBST Warm up RNAscope HiPlex probe stocks and diluent at 40 °C for about 10 min Prepare target probe with warmed HiPlex probe stocks and diluent, mix well and cool to RT before use Set HybEZOven and prepared humidity control tray at 40°C Set RNAscope HiPlex Amp 1-3 and RNAscope HiPlex Fluoro T1-T4 V2 reagents at RT Bring RNAscope HiPlex FFPE Reagent to RT, Prepare 2.5-5% FFPE reagent by using 1:40-1:20 ratio of FFPE reagent to 4X SSC			
14	Probe Incubation	Diluted Probes	2 hrs	40°C
15	Wash X2	Wash Buffer	2 min	RT
16	Amplification	(kit reagents)		
		a) AMP1	30 min	40°C
		b) WASH Buffer	2x2 min	RT
		c) AMP2	30 min	40°C
		d) WASH Buffer	2x2 min	RT
		e) AMP3	30 min	40°C
		f) WASH Buffer	2x2 min	RT
17	FFPE Reagent (optional)	diluted 2.5-5% in 4 X SSC	30 min	RT
18	Wash X2 (optional, only FFPE Reagent)	Wash Buffer	2 min	RT
19	RNAscope HiPlex Fluoro T1-T4 V2	Kit Reagent	15 min	40°C
20	Wash X2	Wash Buffer	2 min	RT
21	Nuclear Staining	DAPI (Kit Reagent)	30 SEC	RT
22	Cover Glass	Prolong Gold		
Store slides in the dark at 2-8 °C				
23	Image slides for round 1	Slide Scanner	6 hrs	
Round 2	Prepare 4 X SSC, 1 L 0.5% PBST Set RNAscope HiPlex Amp 1-3 and RNAscope HiPlex Fluoro T5-T8 V2 reagents and FFPE reagent (Optional) at RT Prepare 10% cleaving solution v1 by diluting with 4X SSC			
24	Remove Coverslips: Soak slides until the cover slips fall off the slides	4 X SSC	1-2 hrs	RT
25	Rinse	4 X SSC	briefly	RT
26	Cleave the T1-T4 V2 fluorophores	10% diluted in 4X SSC	15 min	RT
27	Wash with PBST X2	0.5% PBST	2 min	RT
28	Cleave the T1-T4 V2 fluorophores	10% diluted in 4X SSC	15 min	RT
29	Wash with PBST X2	0.5% PBST	2 min	RT
30	FFPE Reagent (optional)	diluted 2.5-5% in 4 X SSC	30 min	RT
31	Wash X2 (optional, only for FFPE Reagent)	Wash Buffer	2 min	RT
32	RNAscope HiPlex Fluoro T5-T8 V2	Kit Reagent	15 min	40°C
33	Wash X2	Wash Buffer	2 min	RT
34	Nuclear Staining	DAPI (Kit Reagent)	30 SEC	RT
35	Cover Glass	Prolong Gold		
Store slides in the dark at 2-8 °C				
36	Image slides for round 2	Slide Scanner	6 hrs	
Round 3	Set RNAscope HiPlex Amp 1-3 and RNAscope HiPlex Fluoro T9-T12 reagents and FFPE reagent (Optional) at RT Prepare 10% cleaving solution v1 by diluting with 4X SSC			
37	Remove Coverslips: Soak slides until the cover slips fall off the slides	4 X SSC	1-2 hrs	RT
38	Rinse	4 X SSC	briefly	RT
39	Cleave the T5-T8 V2 fluorophores	10% diluted in 4X SSC	15 min	RT
40	Wash with PBST X2	0.5% PBST	2 min	RT
41	Cleave the T5-T8 V2 fluorophores	10% diluted in 4X SSC	15 min	RT
42	Wash with PBST X2	0.5% PBST	2 min	RT
43	FFPE Reagent (optional)	diluted 2.5-5% in 4 X SSC	30 min	RT
44	Wash X2 (optional, only for FFPE Reagent)	Wash Buffer	2 min	RT
45	RNAscope HiPlex Fluoro T9-T12 V2	Kit Reagent	15 min	40°C
46	Wash X2	Wash Buffer	2 min	RT
47	Nuclear Staining	DAPI (Kit Reagent)	30 SEC	RT

Coplín Jars

In Staining Reck

Laying Flat

48	Cover Glass	Prolong Gold	
Store slides in the dark at 2-8 °C			
49	Image slide for round 3	Slide Scanner	6 hrs

Fixed Frozen Tissues:

Fix samples

1. If needed, perfuse tissue with freshly prepared 4% paraformaldehyde (PFA) in 1X PBS or go directly to Step 2.

Note: We recommend perfusing tissues with 1X PBS followed by freshly prepared 4% paraformaldehyde (PFA) in 1X PBS. For suboptimally prepared samples, you may need to adjust pretreatment conditions.

2. Dissect tissue and fix in freshly prepared 10% NBF or 4% PFA for 24 HRS at 4°C.

Freeze tissues

1. Immerse the tissue in 10% sucrose in 1X PBS at 4°C until the tissue sinks to the bottom of the container (approximately 18 HRS for brain tissue).

2. Repeat this step with 20% sucrose in 1X PBS, followed by 30% sucrose in 1X PBS, each time allowing the tissue to sink to the bottom of the container. Freeze the tissue in Optimal Cutting Temperature (OCT) embedding media with dry ice, liquid nitrogen, or 2-methyl butane, and store it in an airtight container at -80°C.

Prepare sections

1. Before tissue sectioning, equilibrate the tissue blocks at -20°C for at least 1 HR in a cryostat.

2. Section the blocks by cutting 7-15 µm thick sections. Mount the sections on SUPERFROST PLUS SLIDES. Do not place sections too close to the edges of the slide.

3. Air dry the slides for 60-120 MIN at -20°C. Use sectioned tissue within three months. Store sections with desiccants at -80°C.