

RNAScope HiPlex Experiment Checklist for FRESH FROZEN Samples				
Method: RNAScope HiPlex				
Experiment Name			Sample Info	
DATE:				
Sample Preparation Steps				
1	Remove the slides from -80°C, immediately Post-fixation	4% PFA/PBS	60 min	RT
2	Rinse X2	PBS	5 min	RT
3	Dehydration			
		a 50% EtOH (50 ml)	5min	RT
		b 70% EtOH (50 ml)	5min	RT
		c 100% EtOH (50 ml)	5min	RT
		d 100% EtOH (50 ml)	5min	RT
		e 100% EtOH (50 ml)	5min	RT
Slides may be stored in 100% Ethanol at -20°C for 1 week				
4	Air Dry on absorbent paper with section face-up	NA	5min	RT
5	Create hydrophobic barrier, let barrier dry completely	NA	5min	RT
Place a humidifying paper in humidity control tray and wet completely with distilled water				
6	Load slides into ACD EZ-Batch Slide Holder	NA	Briefly	RT
7	Protease Digestion	Protease IV (kit)	30min	RT+hum
8	Rinse X2	PBS	2 min	RT
HiPlex Protocol				
Round 1 Prepare 3 L 1 X Wash Buffer				
Warm up RNAScope HiPlex probe stocks and diluent at 40 °C for about 10 min				
Prepare target probe with warmed HiPlex probe stocks and diluent, mix well and cool to RT before use				
Set HybEZoven and prepared humidity control tray at 40°C				
Set RNAScope HiPlex Amp 1-3 and RNAScope HiPlex Fluoro T1-T4 V2 reagents at RT				
Bring RNAScope HiPlex FFPE Reagent to RT, Prepare 2.5-5% FFPE reagent by using 1:40-1:20 ratio of FFPE reagent to 4X SSC				
9	Probe Incubation	Diluted Probes	2 hrs	40°C
10	Wash X2	Wash Buffer	2 min	RT
11	Amplification	(kit reagents)		
		a) AMP1	30 min	40°C
		b) <i>WASH Buffer</i>	2x2 min	RT
		c) AMP2	30 min	40°C
		d) <i>WASH Buffer</i>	2x2 min	RT
		e) AMP3	30 min	40°C
		f) <i>WASH Buffer</i>	2x2 min	RT
12	FFPE Reagent (optional)	diluted 2.5-5% in 4 X SSC	30 min	RT
13	Wash X2 (optional, only FFPE Reagent)	Wash Buffer	2 min	RT
14	RNAscope HiPlex Fluoro T1-T4 V2	Kit Reagent	15 min	40°C
15	Wash X2	Wash Buffer	2 min	RT
16	Nuclear Staining	DAPI (Kit Reagent)	30 SEC	RT
15	Cover Glass	Prolong Gold		
Store slides in the dark at 2-8 °C				
16	Image slides for round 1	Slide Scanner	6 hrs	
Round 2 Prepare 500 ml 4 X SSC, 1 L 0.5% PBST				
Set RNAScope HiPlex Amp 1-3 and RNAscope HiPlex Fluoro T5-T8 V2 reagents and FFPE reagent (Optional) at RT				
Prepare 10% cleaving solution v1 by diluting with 4X SSC				
17	Remove Coverslips: Soak slides until the cover slips fall off the slides	4 X SSC	1-2 hrs	RT
18	Rinse	4 X SSC	briefly	RT
19	Cleave the T1-T4 V2 fluorophores	10% diluted in 4X SSC	15 min	RT
20	Wash with PBST X2	0.5% PBST	2 min	RT
19	Cleave the T1-T4 V2 fluorophores	10% diluted in 4X SSC	15 min	RT
20	Wash with PBST X2	0.5% PBST	2 min	RT
21	FFPE Reagent (optional)	diluted 2.5-5% in 4 X SSC	30 min	RT
22	Wash X2 (optional, only for FFPE Reagent)	Wash Buffer	2 min	RT
23	RNAscope HiPlex Fluoro T5-T8 V2	Kit Reagent	15 min	40°C
24	Wash X2	Wash Buffer	2 min	RT
25	Nuclear Staining	DAPI (Kit Reagent)	30 SEC	RT
26	Cover Glass	Prolong Gold		
Store slides in the dark at 2-8 °C				
27	Image slides for round 2	Slide Scanner	6 hrs	
Round 3 Prepare 1 L 0.5% PBST				
Set RNAScope HiPlex Amp 1-3 and RNAscope HiPlex Fluoro T9-T12 reagents and FFPE reagent (Optional) at RT				
Prepare 10% cleaving solution v1 by diluting with 4X SSC				
28	Remove Coverslips: Soak slides until the cover slips fall off the slides	4 X SSC	1-2 hrs	RT
29	Rinse	4 X SSC	briefly	RT
30	Cleave the T5-T8 V2 fluorophores	10% diluted in 4X SSC	15 min	RT
31	Wash with PBST X2	0.5% PBST	2 min	RT
32	Cleave the T5-T8 V2 fluorophores	10% diluted in 4X SSC	15 min	RT
33	Wash with PBST X2	0.5% PBST	2 min	RT
34	FFPE Reagent (optional)	diluted 2.5-5% in 4 X SSC	30 min	RT
35	Wash X2 (optional, only for FFPE Reagent)	Wash Buffer	2 min	RT
36	RNAscope HiPlex Fluoro T9-T12 V2	Kit Reagent	15 min	40°C
37	Wash X2	Wash Buffer	2 min	RT
38	Nuclear Staining	DAPI (Kit Reagent)	30 SEC	RT
39	Cover Glass	Prolong Gold		
Store slides in the dark at 2-8 °C				
40	Image slide for round 3	Slide Scanner	6 hrs	

Coplin Jars

Prepare fresh-frozen tissue sections

1. Remove tissue and cut to fit into cryomolds.
2. Freeze the specimen on dry ice or in liquid nitrogen, isopentane, or 2-methyl butane within 5 MIN of tissue harvest.
3. Embed the frozen tissue in cryo-embedding medium (OCT) or Tissue Freezing Medium (TFM):
 - a. Add two drops of OCT into a cryomold.
 - b. Place the frozen tissue on the OCT in the correct orientation for cutting.
 - c. Add more OCT to fill the cryomold. Do not allow any air bubbles to form.
 - d. Hold the block with forceps on the surface of the liquid nitrogen or isopentane cooled by dry ice or liquid nitrogen, or place the cryomold on dry ice.
4. Store the frozen block in an air-tight container at -80°C prior to sectioning. Embedded tissue may be stored for up to three months.
5. Section the block:
 - e. Equilibrate block to -20°C in a cryostat ~ 1 HR.
 - f. Cut $10\text{--}20\ \mu\text{m}$ thick sections and mount onto SUPERFROST PLUS SLIDES.
 - g. Dry the sections at 60 -- 120 MIN at -20°C to retain tissue adherence.
6. Store the sections in slide boxes wrapped air-tight with aluminum foil or zip-lock bags at -80°C until use.